

# UTP18 Rabbit pAb

货号: **AYP18432**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:200 - 1:2000
理论分子量	62kDa
实测分子量	60-70kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,Mouse brain,Mouse kidney,Mouse liver
细胞定位	Nucleus,nucleolus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 70-260 of human UTP18 (NP_057085.2).
序列	RRLRQRNRLRLEEDKPAVERCLEELVFGDVENEDALLRRLRGPRVQEHEDSGDSEVENEAKGNFPPQKKPVWVDEEDE DEEMVDMNNRFRKDMMKNAESKLSKDNLKKRLKEEFQHAMGGVPAWAETTKRKTSSDDESEDEDDLLQRTGNFI STSTSLPRGILKMKNQCQHANAERPTVARISSVQFH

## 靶点信息

研究背景	Involved in nucleolar processing of pre-18S ribosomal RNA.
基因ID	51096
基因名	UTP18

<b>Swiss</b>	Q9Y5J1
别名	UTP18;CGI-48;WDR50

## 产品验证



Western blot analysis of UTP18 expressed in Mouse brain,HeLa,Jurkat using UTP18 Rabbit pAb at 1:1000 . Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)