

RBPMS Rabbit pAb

货号: **AYP18462**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	15kDa/17kDa/21kDa/22kDa/24kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart,Mouse kidney,Mouse liver,Mouse lung,Rat testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-219 of human RBPMS (NP_001008712.1).
序列	LVGTPNPSTPLPNTVPQFIAREPYELTVPALYPSSPEVWAPYPLYPALAPALPPPAFTYPASLHAQCFSPPEAKPNTVPFCPLLQQIRFVSGNVFVTYQPTADQQRELPC

靶点信息

研究背景	This gene encodes a member of the RNA recognition motif family of RNA-binding proteins. The RNA recognition motif is between 80-100 amino acids in length and family members contain one to four copies of the motif. The RNA recognition motif consists of two short stretches of conserved sequence, as well as a few highly conserved hydrophobic residues. The encoded protein has a single, putative RNA recognition motif in its N-terminus. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	11030

基因名	RBPM5
Swiss	Q93062
别名	RBPM5;HERMES

产品验证



Western blot analysis of RBPM5 expressed in Mouse heart,Rat testis,Mouse kidney using RBPM5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn