

ZSCAN21 Rabbit pAb

货号: **AYP18501**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	53kDa
实测分子量	54kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,293T,Mouse brain,Mouse heart,Rat brain,Rat heart
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 130-280 of human ZSCAN21 (NP_666019.1).
序列	PGHQVSTPPNEQKPVWEKISSSGTAKESPSSMQPQPLETSHKYESWGPLYIQESGEEQEFAQDPRKVRDCRLSTQHEES ADEQKGSEAEGLKGDIIISVIIANKPEASLERQCVNLENEKGTKPPLQEAGSKKGRESVPTKPTPGERRYICA

靶点信息

研究背景	Strong transcriptional activator (By similarity. Plays an important role in spermatogenesis; essential for the progression of meiotic prophase I in spermatocytes (By similarity).
基因ID	7589
基因名	ZSCAN21

Swiss	Q9Y5A6
别名	ZSCAN21;NY-REN-21;ZNF38;Zipro1

产品验证



Western blot analysis of ZSCAN21 expressed in Mouse brain,Rat brain,MCF7 using ZSCAN21 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120 S.

实验步骤

访问官网浏览详情: www.ablybio.cn