

TNR Rabbit pAb

货号: **AYP18513**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:200 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	139kDa/149kDa
实测分子量	149kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SH-SY5Y,Mouse brain,Rat brain
细胞定位	Secreted,extracellular matrix,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1059-1358 of human TN R (NP_003276.3).
序列	ISWQPPRAEIQENYVLTYSKSTDGSRKELIVDAEDTWIRLEGLLENTDYTVLLQAAQDTTWSSITSTAFTTGGRVFPHPQDCAQHLMNGDTLSGVYPIFLNGELSQKLQVYCDMTDGGGWIVFQRRQNGQTDFFRKWADYRVGFGNVEDEFWLGLDNIHRI TSQGRYELRVDMRDQEAFAFASYDRFSVEDSRNLYKLRIGSYNGTAGDSLSYHQGRPFSTEDRDNDVAVTNCAMSYKG AWWWKNCHRNLNGKYGESRHSQGINWYHWKGHEFSIPFVEMKMRPYNHRLMAGRKRQSLQF

靶点信息

研究背景	This gene encodes a member of the tenascin family of extracellular matrix glycoproteins. The encoded protein is restricted to the central nervous system. The protein may play a role in neurite outgrowth, neural cell adhesion and modulation of sodium channel function. It is a constituent of perineuronal nets.
------	---

基因ID	7143
基因名	TNR
Swiss	Q92752
别名	TNR;TN-R

产品验证



Western blot analysis of TNR expressed in SH-SY5Y, Mouse brain, Rat brain using TNR Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn