

STAC Rabbit pAb

货号: **AYP18515**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:200 - 1:2000
理论分子量	44kDa
实测分子量	44kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,293T,A431,LO2,Mouse brain,Mouse heart,Mouse liver,Rat kidney
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 153-402 of human STAC (NP_003140.1).
序列	GLAPQRCMGKLPKGFRRYYSSPLLIHEQFGCIKEVMPIACGNKVDPVYETLRFGTSLAQRTKKGSSGSGSDSPHRTSTSDL VEVPEEANGPGGGYDLRKRNSVFTYPENGTDDFRDPAKNINHQGSLKDPQMNTYVALYKFVPQENEDLEMRPGDIIT LLEDNEDWWKGIQDRIGFFPANFVQRLQNEKIFRCVTRFIGCKEQGITLKENQICVSSEEEQDGFIRVLSGKKKGLIP LDVLENI

靶点信息

研究背景	Promotes expression of the ion channel CACNA1H at the cell membrane, and thereby contributes to the regulation of channel activity. Plays a minor and redundant role in promoting the expression of calcium channel CACNA1S at the cell membrane, and thereby contributes to increased channel activity. Slows down the inactivation rate of the calcium channel CACNA1C.
------	---

基因ID	6769
基因名	STAC
Swiss	Q99469
别名	STAC;STAC1

产品验证



Western blot analysis of STAC expressed in Mouse brain, Rat kidney, U-251MG using STAC Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn