

— ABLYBIO, Help Your Research



# CEBP Delta/CEBPD Rabbit pAb

货号: **AYP18552**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	28kDa
实测分子量	28KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	3T3L1,RAW264.7,Mouse brain,Rat brain
细胞定位	Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-150 of human CEBP Delta/CEBPD (NP_005186.2).
------	--

## 靶点信息

研究背景	The protein encoded by this intronless gene is a bZIP transcription factor which can bind as a homodimer to certain DNA regulatory regions. It can also form heterodimers with the related protein CEBP-alpha. The encoded protein is important in the regulation of genes involved in immune and inflammatory responses, and may be involved in the regulation of genes associated with activation and/or differentiation of macro phages. The cytogenetic location of this locus has been reported as both 8p11 and 8q11.
基因ID	1052
基因名	CEBPD
Swiss	P49716 ( <a href="https://www.uniprot.org/uniprotkb/P49716/entry">https://www.uniprot.org/uniprotkb/P49716/entry</a> )
别名	CEBPD,C/EBP-delta,CELF,CRP3,NF-IL6-beta,CEBP Delta/CEBPD Rabbit pAb,Nuclear factor NF-IL6-beta

## 产品验证



Western blot analysis of CEBP Delta/CEBPD expressed in Mouse brain,Rat brain,3T3L1 using CEBP Delta/C EBPD Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/prote ins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure ti me: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)