

ANGPTL5 Rabbit pAb

货号: **AYP18576**

产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	44kDa
实测分子量	44kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Rat kidney, Rat heart
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-388 of human ANGPTL5 (NP_835228.2).
序列	MMSPSQASLLFLNVCIFICGEAVQGNVCVHHSTDSVVNIVEDGSSNAKDESKSNDTVCKEDCEESCDVTKITREEKHFMC RNLQNSIVSYTRSTKLLRNMMDEQQASLDYLSNQVNELMNRVLLLTTEVFRKQLDPFPHRPVQSHGLDCTDIKDTIGSVT KTPSGLYIIHPEGSSYPFEVMCDMDYRGGGRTVIQKRIDGIIDFQRLWCDYLDGFGDLLGEFWLGLKKIFYIVNQKNTSFML YVALESEDDTLAYASYDNFWLEDETRFFKMHLGRYSGNAGDAFRGLKKEDNQNAMPFSTSDVDNDGCRPACLVNGQS VKSCSHLHNKTGWVWFNECGLANLNGIHHFSGKLLATGIQWGTWTKNNSPVKIKSVSMKIRRMYPYFK

靶点信息

研究背景	Predicted to be active in collagen-containing extracellular matrix and extracellular space.
基因ID	253935

基因名	ANGPTL5
Swiss	Q86XS5
别名	ANGPTL5

产品验证



Western blot analysis of ANGPTL5 expressed in Rat kidney,Rat heart using ANGPTL5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn