

— ABLYBIO, Help Your Research



# SLC44A2 Rabbit pAb

货号: **AYP18596**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	79kDa/80kDa
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SGC-7901,Mouse lung,Mouse thymus,Rat lung
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 150-250 of human SLC44A2 (NP_065161.3).
------	--

## 靶点信息

研究背景	Isoform 1, but not isoform 3, exhibits some choline transporter activity.
------	---

基因ID	57153
基因名	SLC44A2
Swiss	Q8IWA5
别名	SLC44A2,CTL2,PP1292,SLC44A2 Rabbit pAb,Solute carrier family 44 member 2

## 产品验证



Western blot analysis of SLC44A2 expressed in Mouse lung,Rat lung,SGC-7901 using SLC44A2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)