

REC8 Rabbit pAb

货号: **AYP18679**

产品信息

反应	Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	67kDa
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Rat testis
细胞定位	Chromosome, Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 194-561 of mouse REC8 (NP_064386.2).
序列	EVITLQEAEPIRMQLQIEGEQDLPEISRGDLELLIAEKDDAILLEERQGRLLRQRASLPLDESREEPRALEGAGLVLSALSPAP AQVEGIQEALPGQVFPPEVQKMTGWEPGALLTEVTPPQELRLPAPPSTEKRLPSLQRPLPRRHRRRQLLFWDKETQISREK FEEQLQTGAHCWEYPVAQPPKRMILTSPAELFRTPTLSGWLPELLGLWTHCAQVPQRMLRQRPQLETEETVEEERAADDEE ERRKTEALSEIEVLREAQEPSGPLMLSSELSLEAAEDEKSRSLIPPEWWAWSEEQPEPPALPMLPELPEVPMEMPPRPEL SSEAVLRAVALKLQANKELDFSSLPPLSPRKLASRVFYL

靶点信息

研究背景	This gene encodes a member of the kleisin family of SMC (structural maintenance of chromosome) protein partners. The protein localizes to the axial elements of chromosomes during meiosis in both oocytes and spermatocytes. In the mouse, the homologous protein is a key component of the meiotic cohesion complex, which regulates sister chromatid cohesion and recombination between homologous chromosomes. Multiple alternatively spliced variants, encoding the same protein, have been found for this gene.
------	---

基因ID	56739
基因名	
Swiss	
别名	REC8;HR21spB;REC8L1;Rec8p

产品验证



Western blot analysis of REC8 expressed in Rat testis using REC8 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn