

— ABLYBIO, Help Your Research



ADAM32 Rabbit pAb

货号: **AYP18693**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 87kDa |
| 实测分子量 | 88kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa,293T,Mouse brain |
| 细胞定位 | Membrane,Single-pass type I membrane protein |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|---|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 703-787 of human ADA M32 (NP_659441.3). |
|------|---|

靶点信息

| | |
|-------|--|
| 研究背景 | This gene encodes a member of the disintegrin family of membrane-anchored proteins that play a role in diverse biological processes such as brain development, fertilization, tumor development and inflammation. This gene is predominantly expressed in the testis. The encoded protein undergoes proteolytic processing to generate a mature polypeptide comprised of an metalloprotease, disintegrin and epidermal growth factor-like domains. This gene is located in a cluster of other disintegrin and metallopeptidase family genes on chromosome 8. Alternative splicing results in multiple transcript variants encoding different isoforms. |
| 基因ID | 203102 |
| 基因名 | ADAM32 |
| Swiss | Q8TC27 (https://www.uniprot.org/uniprotkb/Q8TC27/entry) |
| 别名 | ADAM32,ADAM32 Rabbit pAb |

产品验证



Western blot analysis of ADAM32 expressed in HeLa,293T,Mouse brain using ADAM32 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)