

LENG8 Rabbit pAb

货号: **AYP18706**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	5kDa/86kDa/88kDa
实测分子量	86kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Jurkat,Mouse brain,Mouse eye,Rat kidney,Rat brain
细胞定位	nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 101-319 of human LENG 8 (NP_443157.1).
序列	YYYPMSMYQSYGSPSQYGMAGSYGSATPQQPSAPQHQTGNQPPVPGMDESMYSYQAPPQQLPSAQPQPSPNPPHGAH TLNSGPQPGTAPATQHSQAGPATGQAYGPHTYTEPAKPKKGQQLWNRMKPAPGTGGLKFNIQKRPFAVTTQSFGSNAE GQHSFGPQPNEKVNHSYGSSARGNLSGKPDWQDMKEYVERCFTACESEEDKDRTEKLLKE

靶点信息

研究背景	Predicted to be active in nucleus.
基因ID	114823
基因名	LENG8

Swiss	Q96PV6
别名	LENG8;pp13842

产品验证



Western blot analysis of LENG8 expressed in Mouse brain,Rat kidney,U-251MG using LENG8 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn