

— ABLYBIO, Help Your Research



DEXI Rabbit pAb

货号: **AYP18748**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	10kDa
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	A-549, BxPC-3, HT-29
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-95 of human DEXI (NP_054734.2).
------	---

靶点信息

研究背景	DEXI (Dexi Homolog) is a Protein Coding gene. Diseases associated with DEXI include Fibrosclerosis Of Breast and Immunoglobulin A Deficiency 1.
------	---

基因ID	28955
基因名	DEXI
Swiss	O95424 (https://www.uniprot.org/uniprotkb/O95424/entry)
别名	DEXI,MYLE,DEXI Rabbit pAb,Protein MYLE

产品验证



Western blot analysis of DEXI expressed in A-549, BxPC-3, HT-29 using DEXI Rabbit pAb at 1:1000. Second ary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)