

— ABLYBIO, Help Your Research



# ITGB3BP Rabbit pAb

货号: **AYP18762**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	IHC
推荐浓度	<b>IHC:</b> 1:50 - 1:200
理论分子量	6kDa/12kDa/19kDa/20kDa/24kDa
实测分子量	6kDa/12kDa/19kDa/20kDa/24kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	Chromosome,Cytoplasm,Cytoplasm,Nucleus,centromere,kinetochore
纯化	Affinity purification

## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-111 of human ITGB3BP (NP_055103.3).
------	---

## 靶点信息

研究背景	This gene encodes a transcriptional coregulator that binds to and enhances the activity of members of the nuclear receptor families, thyroid hormone receptors and retinoid X receptors. This protein also acts as a corepressor of NF-kappaB-dependent signaling. This protein induces apoptosis in breast cancer cells through a caspase 2-mediated signaling pathway. This protein is also a component of the centromere-specific histone H3 variant nucleosome associated complex (CENP-NAC) and may be involved in mitotic progression by recruiting the histone H3 variant CENP-A to the centromere. Alternate splicing results in multiple transcript variants.
基因ID	23421
基因名	ITGB3BP
Swiss	Q13352 ( <a href="https://www.uniprot.org/uniprotkb/Q13352/entry">https://www.uniprot.org/uniprotkb/Q13352/entry</a> )
别名	ITGB3BP,CENP-R,CENPR,HSU37139,NRIF3,TAP20,ITGB3BP Rabbit pAb,Beta-3-endonexin,Integrin beta-3-binding protein,Nuclear receptor-interacting factor 3

## 产品验证

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)