

— ABLYBIO, Help Your Research



ARL6IP5 Rabbit pAb

货号: **AYP18768**

产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	21kDa
实测分子量	22kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse thymus
细胞定位	Cytoplasm, Endoplasmic reticulum membrane, Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100 to the C-terminus of human ARL6IP5 (NP_006398.1).
------	--

靶点信息

研究背景	Expression of this gene is affected by vitamin A. The encoded protein of this gene may be associated with the cytoskeleton. A similar protein in rats may play a role in the regulation of cell differentiation. The rat protein binds and inhibits the cell membrane glutamate transporter EAAC1. The expression of the rat gene is upregulated by retinoic acid, which results in a specific reduction in EAAC1-mediated glutamate transport.
基因ID	10550
基因名	ARL6IP5
Swiss	O75915 (https://www.uniprot.org/uniprotkb/O75915/entry)
别名	ARL6IP5,DERP11,GTRAP3-18,HSPC127,JWA,PRAF3,Yip6b,addicisin,hp22,jmx,ARL6IP5 Rabbit pAb,ADP-ribosylation factor-like protein 6-interacting protein 5,Cytoskeleton-related vitamin A-responsive protein,Dermal papilla-derived protein 11,Glutamate transporter EAAC1-interacting protein,JM5

产品验证



Western blot analysis of ARL6IP5 expressed in Mouse thymus using ARL6IP5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)