

FAM127A Rabbit pAb

货号: **AYP18781**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	
实测分子量	15kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cell membrane,Lipid-anchor
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-113 of human FAM127A (NP_001071639.1).
序列	MDGRVQLIKALLALPIRPATRRWRNPIPFPEFDGDTDRLPEFIVQTGSYMFVDENTFSSDALKVTFILITRLTGPALQWVIPYIKKESPLLNDYRGFLAEMKRVFGWEEDEDF

靶点信息

研究背景	Predicted to be located in plasma membrane.
基因ID	8933
基因名	RTL8C

Swiss	A6ZK13
别名	RTL8C;CXX1;FAM127A;MAR8C;MART8C;Mar8;Mart8;SIRH5

产品验证



Western blot analysis of FAM127A expressed in Mouse brain,Rat brain using FAM127A Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn