

PRPF4B Rabbit pAb

货号: **AYP18784**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	116kDa
实测分子量	160kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,A-549,HeLa,Mouse liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 550-760 of human PRPF 4B (NP_003904.3).
序列	QKYKYLAEDSNMSVPSEPSSPQSSTRTRSPSPDDILERVAADVKEYERENVDTFEASVKAKHNLMTVEQNGSSQKLLA PDMFTESDDMFAAYFDSARLRAAGIGKDFKENPNLRDNWTD AEGYYRVNIGEVLDKRYNVVGYTGQGVFSNVVRARDN ARANQEVAVKIIRNNELMQKTGLKELEFLKLLNDADPDDKFHCLRLFRHFYHK

靶点信息

研究背景	Pre-mRNA splicing occurs in two sequential transesterification steps, and the protein encoded by this gene is thought to be involved in pre-mRNA splicing and in signal transduction. This protein belongs to a kinase family that includes serine/arginine-rich protein-specific kinases and cyclin-dependent kinases (CDKs). This protein is regarded as a CDK-like kinase (Clk) with homology to mitogen-activated protein kinases (MAPKs).
------	--

基因ID	8899
基因名	PRPF4B
Swiss	Q13523
别名	PRPF4B;PR4H;PRP4;PRP4H;PRP4K;dj1013A10.1

产品验证



Western blot analysis of PRPF4B expressed in Mouse liver,293T,A-549 using PRPF4B Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn