

FOXS1 Rabbit pAb

货号: **AYP18825**

产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:100 - 1:500
理论分子量	35kDa
实测分子量	33kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Rat testis
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 121-330 of human FOXS 1 (NP_004109.1).
序列	EGTRGPAKARRGPLRATSQDPGVNATTGRQCSFPPELPDPKGLSFGGLVGMAMPASMCPATTDGRPRPPMEPKEISTPKP ACPGELPVATSSSSCPAFGFPAGFSEAESFNKAPTPLVSPESGIGSSYQCRLQALNFCMGADPGLEHLLASAAPSPAPPTPP GSLRAPLPLPTDHKEPWVAGGFPVQGGSGYPLGLTPCLYRTPGMFFFE

靶点信息

研究背景	The forkhead family of transcription factors belongs to the winged helix class of DNA-binding proteins. The e protein encoded by this intronless gene contains a forkhead domain and is found predominantly in aort a and kidney. The function of the encoded protein is unknown.
基因ID	2307

基因名	FOXS1
Swiss	O43638
别名	FOXS1;FKHL18;FREAC10

产品验证



Western blot analysis of FOXS1 expressed in 293T,Rat testis using FOXS1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn