

# HIF1 $\beta$ /ARNT Rabbit pAb

货号: **AYP18844**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:1000 - 1:5000
理论分子量	84kDa/85kDa/86kDa
实测分子量	87KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat uterus
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 465-570 of human HIF1 $\beta$ /ARNT (NP_001659.1).
序列	KNSSQEPRPTLSNTIQRPQLGPTANLPLEMGSQLAPRQQQQTELDMPGRDGLASYNHSQVVQPVTITGPEHSKPLE KSDGLFAQDRDPRFSEIYHNINADQSK

## 靶点信息

研究背景	This gene encodes a protein containing a basic helix-loop-helix domain and two characteristic PAS domains along with a PAC domain. The encoded protein binds to ligand-bound aryl hydrocarbon receptor and aids in the movement of this complex to the nucleus, where it promotes the expression of genes involved in xenobiotic metabolism. This protein is also a co-factor for transcriptional regulation by hypoxia-inducible factor 1. Chromosomal translocation of this locus with the ETV6 (ets variant 6) gene on chromosome 12 have been described in leukemias. Alternative splicing results in multiple transcript variants.
------	---

基因ID	405
基因名	ARNT
Swiss	P27540
别名	ARNT;HIF-1-beta;HIF-1beta;HIF1-beta;HIF1B;HIF1BETA;TANGO;bHLHe2

## 产品验证



Western blot analysis of HIF1 $\beta$ /ARNT expressed in Rat uterus using HIF1 $\beta$ /ARNT Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)