

# IGF1R Rabbit pAb

货号: **AYP18846**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	154kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Mouse testis,Rat testis
细胞定位	Cell membrane,Single-pass type I membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 741-910 of human IGF1 R (NP_000866.1).
序列	DVMQVANTTMSRSRNTTAADTYNITDPEELETEYPPFESRVDNKERTVISNLRPFTLYRIDIHSCNHEAEKLGCSASNFVF ARTMPAEGADDIPGPVTWEPRPENSIFLKWPEPENGLILMYEIKYGSQVEDQRECVSRQEYRKYGGAKLNRLNPGNYTA RIQATSL

## 靶点信息

研究背景	This receptor binds insulin-like growth factor with a high affinity. It has tyrosine kinase activity. The insulin-like growth factor I receptor plays a critical role in transformation events. Cleavage of the precursor generates alpha and beta subunits. It is highly overexpressed in most malignant tissues where it functions as an anti-apoptotic agent by enhancing cell survival. Alternatively spliced transcript variants encoding distinct isoforms have been found for this gene.
------	---

基因ID	3480
基因名	IGF1R
Swiss	P08069
别名	IGF1R;CD221;IGFIR;IGFR;JTK13

## 产品验证



Western blot analysis of IGF1R expressed in U-251MG, Mouse testis, Rat testis using IGF1R Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)