

— ABLYBIO, Help Your Research



DCSTAMP Rabbit pAb

货号: **AYP18869**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:5000
理论分子量	31kDa/53kDa
实测分子量	53KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain
细胞定位	Cell membrane,Endoplasmic reticulum membrane,Endoplasmic reticulum-Golgi intermediate compartment membrane,Endosome,Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 401-470 of human DCSTAMP (NP_110415.1).
------	---

靶点信息

研究背景	This gene encodes a seven-pass transmembrane protein that is primarily expressed in dendritic cells. The encoded protein is involved in a range of immunological functions carried out by dendritic cells. This protein plays a role in osteoclastogenesis and myeloid differentiation. Alternate splicing results in multiple transcript variants.
基因ID	81501
基因名	DCSTAMP
Swiss	Q9H295 (https://www.uniprot.org/uniprotkb/Q9H295/entry)
别名	DCSTAMP, FIND, TM7SF4, hDC-STAMP, DCSTAMP Rabbit pAb, Dendrocyte-expressed seven transmembrane protein, IL-four-induced protein, Transmembrane 7 superfamily member 4

产品验证



Western blot analysis of DCSTAMP expressed in Mouse brain using DCSTAMP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)