

# KLHL1 Rabbit pAb

货号: **AYP18887**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	82kDa
实测分子量	83kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human KLHL1 ( NP_065917.1).
序列	MSGSGRKDFDVKHILRLRWKLFSHPSPTGGPAGGGCLQQDGGSGSFEHWGSPQSRLKSKQERSGVSTFWKKPSSSSSS SSSPSSSSSSFNPLNGTLLPVATRLQQGAPGQGTQQPARTLFYVESLEEEVVPGMDFPGPHEKGLVLQELKVEPDNSSQA TGEGCGHRLSSTGHSMTQPQSDLSSSSEEFYQAVHHAETFRKMESYLKQQQ

## 靶点信息

研究背景	
基因ID	57626
基因名	KLHL1

Swiss	Q9NR64
别名	KLHL1;MRP2

## 产品验证



Western blot analysis of KLHL1 expressed in A-431 using KLHL1 Rabbit pAb at 1:1000. Secondary antibody y: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)