

— ABLYBIO, Help Your Research



# SNRPD1 Rabbit pAb

货号: **AYP18933**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	13kDa
实测分子量	13KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,K-562,Mouse brain,Mouse spleen,Rat brain,Rat thymus
细胞定位	Cytoplasm,Nucleus,cytosol
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-96 of human SNRPD1 (NP_008869.1).
------	---

## 靶点信息

研究背景	This gene encodes a small nuclear ribonucleoprotein that belongs to the SNRNP core protein family. The protein may act as a charged protein scaffold to promote SNRNP assembly or strengthen SNRNP-SNRNP interactions through nonspecific electrostatic contacts with RNA. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	6632
基因名	SNRPD1
Swiss	P62314 ( <a href="https://www.uniprot.org/uniprotkb/P62314/entry">https://www.uniprot.org/uniprotkb/P62314/entry</a> )
别名	SNRPD1,HsT2456,SMD1,SNRPD,Sm-D1,SNRPD1 Rabbit pAb,Sm-D autoantigen,snRNP core protein D1

## 产品验证



Western blot analysis of SNRPD1 expressed in Mouse brain,Rat brain,HeLa using SNRPD1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))