

— ABLYBIO, Help Your Research



FAM117A Rabbit pAb

货号: **AYP18936**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	48.3kDa
实测分子量	48KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,K-562,Mouse liver
细胞定位	
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-240 of human FAM117A (NP_110429.1).
------	---

靶点信息

研究背景	FAM117A (Family With Sequence Similarity 117 Member A) is a Protein Coding gene. An important paralog of this gene is FAM117B.
基因ID	81558
基因名	FAM117A
Swiss	Q9C073 (https://www.uniprot.org/uniprotkb/Q9C073/entry)
别名	FAM117A Rabbit pAb,FAM117A,C/EBP-induced protein

产品验证



Western blot analysis of FAM117A expressed in HepG2, K-562, Mouse liver using FAM117A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)