

# FAM117A Rabbit pAb

货号: **AYP18936**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200
理论分子量	48.3kDa
实测分子量	48KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,K-562,Mouse liver
细胞定位	
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 110-240 of human FAM 117A (NP_110429.1).
序列	GAFTCCTNDKATQTPLSWQELEGERASSCAHKRSASWGSTDHRKEISKLKQQLQRTKLSRSGKEKERGSPLLGDHAVR GALRASPPSPSGSPVLRSLSPCLHRSLEGLNQELEEVFVKEQGEEELLRILDI

## 靶点信息

研究背景	FAM117A (Family With Sequence Similarity 117 Member A) is a Protein Coding gene. An important paralog of this gene is FAM117B.
基因ID	81558

基因名	FAM117A
Swiss	Q9C073
别名	

## 产品验证



Western blot analysis of FAM117A expressed in HepG2, K-562, Mouse liver using FAM117A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)