

— ABLYBIO, Help Your Research



# PEBP4 Rabbit pAb

货号: **AYP18939**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	25kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse lung
细胞定位	Lysosome
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 23-120 of human PEBP4 (NP_659399.2).
------	--

## 靶点信息

研究背景	
------	--

基因ID	157310
基因名	PEBP4
Swiss	Q96S96 ( <a href="https://www.uniprot.org/uniprotkb/Q96S96/entry">https://www.uniprot.org/uniprotkb/Q96S96/entry</a> )
别名	PEBP4,CORK-1,CORK1,GWTM1933,HEL-S-300,PEBP-4,PRO4408,hPEBP4,PEBP4 Rabbit pAb,Protein cousin-of-RKIP 1

## 产品验证



Western blot analysis of PEBP4 expressed in Mouse lung using PEBP4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)