

RXFP2 Rabbit pAb

货号: **AYP18954**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 83kDa/86kDa |
| 实测分子量 | 86kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | U-87MG,293T,HeLa,K-562,Mouse brain,Mouse kidney,Mouse testis,Rat brain,Rat kidney |
| 细胞定位 | Cell membrane,Multi-pass membrane protein |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human RXFP2 (NP_570718.1). |
| 序列 | CGDTSGWATIFGTVHGNANSVALTQECFLKQYPQCCDCKETELEC VNGDLKSVPMISNNVTLLSLKKNKIHS L PDKVFIKY TKLKKIFLQHNCIRHISRKA |

靶点信息

| | |
|------|--|
| 研究背景 | This gene encodes a member of the GPCR (G protein-coupled, 7-transmembrane receptor) family. Mutations in this gene are associated with cryptorchidism. Alternatively spliced transcript variants encoding different isoforms have been found for this gene. |
| 基因ID | 122042 |

| | |
|-------|---|
| 基因名 | RXFP2 |
| Swiss | Q8WXD0 |
| 别名 | RXFP2;GPR106;GREAT;INSL3R;LGR8;LGR8.1;RXFP2 |

产品验证



Western blot analysis of RXFP2 expressed in Mouse brain,Rat brain,U-87MG using RXFP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn