

— ABLYBIO, Help Your Research



RXFP2 Rabbit pAb

货号: **AYP18954**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	83kDa/86kDa
实测分子量	86kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,293T,HeLa,K-562,Mouse brain,Mouse kidney,Mouse testis,Rat brain,Rat kidney
细胞定位	Cell membrane,Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 80-180 of human RXFP2 (NP_570718.1).
------	--

靶点信息

研究背景	This gene encodes a member of the GPCR (G protein-coupled, 7-transmembrane receptor) family. Mutations in this gene are associated with cryptorchidism. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
基因ID	122042
基因名	RXFP2
Swiss	Q8WXD0 (https://www.uniprot.org/uniprotkb/Q8WXD0/entry)
别名	RXFP2,GPR106,GREAT,INSL3R,LGR8,LGR8.1,RXFPR2,RXFP2 Rabbit pAb,G-protein coupled receptor 106,G-protein coupled receptor affecting testicular descent,Leucine-rich repeat-containing G-protein coupled receptor 8,Relaxin family peptide receptor 2

产品验证



Western blot analysis of RXFP2 expressed in Mouse brain,Rat brain,U-87MG using RXFP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)