

# PCNP Rabbit pAb

货号: **AYP18962**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	6kDa/13kDa/18kDa
实测分子量	19kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG, HT-29, HeLa
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-178 of human PCNP (NP_065090.1).
序列	MADGKAGDEKPEKSQRAGAAGGPPEEEAEKPVKTKTVSSSNGGESSRSAEKRSAEAAAADLPTKPTKISKFGFAIGSQT KKASAIKLGSSKPKETVPTLAPKTLVAAAFNEDEDESEPEEMPEAKMRMKNIGRDTPTSAGPNSFNKKGKHGFSQKL WERNIKSHLGNVHDQDN

## 靶点信息

研究背景	May be involved in cell cycle regulation.
基因ID	57092
基因名	PCNP

Swiss	Q8WW12
别名	PCNP

## 产品验证



Western blot analysis of PCNP expressed in U-87MG, HT-29, HeLa using PCNP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)