

RNF181 Rabbit pAb

货号: **AYP18963**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	17kDa
实测分子量	18kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	LO2,Mouse liver,Mouse brain,Mouse kidney,Mouse lung,Rat testis,Rat brain
细胞定位	cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-153 of human RNF181 (NP_057578.1).
序列	MASYFDEHDCEPSDPEQETRNTNMLELARSLFNRMDFEDLGLVVDWDHHLPPPAAKTVVENLPRTVIRGSQAEKCPVCL LEFEEEEETAIEMPCHHLFHSSCILPWLSKTNNSCPLCRYELPTDDDDTYEEHRRDKARKQQQHRLENLHGAMYT

靶点信息

研究背景	RNF181 binds the integrin alpha-IIb (ITGA2B; MIM 607759)/beta-3 (ITGB3; MIM 173470) complex and has E3 ubiquitin ligase activity (Brophy et al., 2008 [PubMed 18331836]).
基因ID	51255
基因名	RNF181

Swiss	Q9P0P0
别名	RNF181;HSPC238

产品验证



Western blot analysis of RNF181 expressed in Mouse liver,Rat testis,LO2 using RNF181 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn