

— ABLYBIO, Help Your Research



# NECAP2 Rabbit pAb

货号: **AYP18986**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	19kDa/25kDa/28kDa/29kDa
实测分子量	28kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney,Mouse liver
细胞定位	Cell membrane,Cytoplasmic vesicle,clathrin-coated vesicle membrane
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 194-263 of human NECAP2 (NP_060560.1).
------	--

## 靶点信息

研究背景	This gene likely encodes a member of the adaptin-ear-binding coat-associated protein family. Studies of a similar protein in rat suggest a role in clathrin-mediated endocytosis. Alternatively spliced transcripts have been described.
基因ID	55707
基因名	NECAP2
Swiss	Q9NVZ3 ( <a href="https://www.uniprot.org/uniprotkb/Q9NVZ3/entry">https://www.uniprot.org/uniprotkb/Q9NVZ3/entry</a> )
别名	NECAP2,NECAP2 Rabbit pAb,NECAP endocytosis-associated protein 2

## 产品验证



Western blot analysis of NECAP2 expressed in Mouse kidney, Mouse liver using NECAP2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)