

# LHPP Rabbit pAb

货号: **AYP18988**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	22kDa/29kDa
实测分子量	29kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse liver,Rat liver,Rat kidney
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-270 of human LHPP (N P_071409.3).
序列	MAPWGKRLAGVRGBVLLDISGVLYDSGAGGGTAIAGSVEAVARLKRSLKVRVFCTNESQKSRAELVGQLQRLGFDISEQEV TAPAPAACQILKEQGLRPYLLIHDGVRSEFDQIDTSNPNCVVIADAGESFSYQNMNNAFQVLMELEKPVLIISLGKGRYYKET SGLMLDVGYPYMKALEYACGIKAEVVGKPSPEFFKSALQAIGVEAHQAVMIGDDIVGDVGGARCGMRALQVRTGKFRPS DEHHPEVKADGYVDNLAEAVDLLLQHADK

## 靶点信息

研究背景	Phosphatase that hydrolyzes imidodiphosphate, 3-phosphohistidine and 6-phospholysine. Has broad substrate specificity and can also hydrolyze inorganic diphosphate, but with lower efficiency (By similarity).
基因ID	64077

基因名	LHPP
Swiss	Q9H008
别名	LHPP;HDHD2B

## 产品验证



Western blot analysis of LHPP expressed in Mouse brain,Rat liver,Mouse liver using LHPP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)