

# TSPAN31 Rabbit pAb

货号: **AYP19028**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	23kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Mouse heart,Mouse kidney,Mouse liver,Rat ovary
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 95-175 of human TSPAN 31 (NP_005972.1).
序列	SCLAINRSKQTDVINASWWVMSNKTRDELERSFDCCGLFNLTLYQQDYDFCTAICKSQSPTCQMCGEKFLKHSDEALKI L

## 靶点信息

研究背景	The protein encoded by this gene is a member of the transmembrane 4 superfamily, also known as the te traspanin family. Most of these members are cell-surface proteins that are characterized by the presence of four hydrophobic domains. The proteins mediate signal transduction events that play a role in the regulation of cell development, activation, growth and motility. This encoded protein is thought to be involved in growth-related cellular processes. This gene is associated with tumorigenesis and osteosarcoma.
基因ID	6302

基因名	TSPAN31
Swiss	Q12999
别名	TSPAN31;SAS

## 产品验证



Western blot analysis of TSPAN31 expressed in Mouse brain,Rat ovary,Mouse heart using TSPAN31 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)