

SERPINB5 Rabbit pAb

货号: **AYP19092**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	25kDa/42kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,NIH/3T3
细胞定位	Secreted,extracellular space
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 21-189 of human SERPINB5 (NP_002630.2).
序列	EKEPLGNVLFSPICLSTLSLAQVGAKGDTANEIGQVLHFENVKDVPPFGFQTVTSDVNKLSSFYSLKLIKRLYVDKSLNLSTEFISSTKRPYAKELETVDKFKLEETKQINNSIKDLTDGHFENILADNSVNDQTKILVVNAAAYFVGKWMKKFSESETKPCFRVVK

靶点信息

研究背景	Tumor suppressor. It blocks the growth, invasion, and metastatic properties of mammary tumors. As it does not undergo the S (stressed to R (relaxed conformational transition characteristic of active serpins, it exhibits no serine protease inhibitory activity.
基因ID	5268

基因名	SERPINB5
Swiss	P36952
别名	SERPINB5;PI5;maspin;serpin B5

产品验证



Western blot analysis of SERPINB5 expressed in HT-29,NIH/3T3 using SERPINB5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn