

— ABLYBIO, Help Your Research



FGL2 Rabbit pAb

货号: **AYP19200**

产品信息

| | |
|-------|------------------------------------------------------------------------------------------------------|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 50kDa |
| 实测分子量 | 70kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HepG2 |
| 细胞定位 | Secreted |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|-------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 24-140 of human FGL2 (NP_006673.1). |
|------|-------------------------------------------------------------------------------------------------------------------|

靶点信息

| | |
|-------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | The protein encoded by this gene is a secreted protein that is similar to the beta- and gamma-chains of fibrinogen. The carboxyl-terminus of the encoded protein consists of the fibrinogen-related domains (FRED). The encoded protein forms a tetrameric complex which is stabilized by interchain disulfide bonds. This protein may play a role in physiologic functions at mucosal sites. |
| 基因ID | 10875 |
| 基因名 | FGL2 |
| Swiss | Q14314 (https://www.uniprot.org/uniprotkb/Q14314/entry) |
| 别名 | FGL2,T49,pT49,FGL2 Rabbit pAb,Fibrinogen-like protein 2 |

产品验证



Western blot analysis of FGL2 expressed in HepG2 using FGL2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)