

DENR Rabbit pAb

货号: **AYP19226**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	22kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,BxPC-3,HeLa,Jurkat,Mouse heart,Rat lung
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-198 of human DENR (NP_003668.2).
序列	MAADISESSGADCKGDPNSAKLDADYPLRVLYCGVCSLPTEYCEYMPDVAKCRQWLEKNFPNEFAKLTVENSPKQEAGI SEGQGTAGEEEEEKKKQKRGGRGQIKQKKKTVPQKVTIKIPRAKKKYVTRVCGLATFEIDLKEAQRFFAQKFSCGASVTGE DEIIIQGDFTDDIIDVVIQEKWPEVDDDSIEDLGEVKK

靶点信息

研究背景	This gene encodes a protein whose expression was found to increase in cultured cells at high density but not during growth arrest. This gene was also shown to have increased expression in cells overexpressing HER-2/neu proto-oncogene. The protein contains an SUI1 domain. In budding yeast, SUI1 is a translation initiation factor that along with eIF-2 and the initiator tRNA-Met, directs the ribosome to the proper translation start site. Proteins similar to SUI have been found in mammals, insects, and plants.
------	---

基因ID	8562
基因名	DENR
Swiss	O43583
别名	DENR;DRP;DRP1;SMAP-3

产品验证



Western blot analysis of DENR expressed in Mouse heart,Rat lung,U-87MG using DENR Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn