

PRDM14 Rabbit pAb

货号: **AYP19248**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	64kDa
实测分子量	70KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NCCIT,A-549
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human PRDM14 (NP_078780.1).
序列	MALPRPSEAVPQDKVCYPPESSPQNLAAYYTPFPSYGHYRNSLATVEEDFQPFQRQLEAAASAAPAMPPFPFRMAPLLSPGLGLQREPLYDLPWYSKLPWYPIPHVPREVPPFLSSSHEYAGASSEDLGHQIIGGDNESGPCGPDTLIPPPADASLLPEGLRTSQQLLPCSPSKQSEDGPKPSNQEGKSPARFQFTEEDLHFVLYGVTPSLEHPASLHHAISGLLVPPDSSGSDSLPQTLDKDSLQLEGLCLMQTVF

靶点信息

研究背景	This gene encodes a member of the PRDI-BF1 and RIZ homology domain containing (PRDM) family of transcriptional regulators. The encoded protein may possess histone methyltransferase activity and plays a critical role in cell pluripotency by suppressing the expression of differentiation marker genes. Expression of this gene may play a role in breast cancer.
------	---

基因ID	63978
基因名	PRDM14
Swiss	Q9GZV8
别名	PRDM14;PFM11

产品验证



Western blot analysis of PRDM14 expressed in NCCIT,A-549 using PRDM14 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn