

FIBP Rabbit pAb

货号: **AYP19280**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	41kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,Mouse skeletal muscle
细胞定位	Endomembrane system,Nucleus,Peripheral membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-364 of human FIBP (N P_004205.2).
序列	MTSELDIFVGNNTLLIDEDVYRLWLDGYSVTDVAVALRVRSGILEQTGATAAVLQSDTMDHYRTFHMLERLLHAPPKLLHQLIF QIPPSRQALLIERYAFDEAFVREVLGKKLSKGTKKDLDDISTKTGITLKCRRQFDNFKRVFKVVEEMRGS�VDNIQQHFL LSDRLARDYAAIVFFANNRFETGKKLQYLSFGDFAFCAELMIQNWTLGAVDSQMDDMDMDLDKEFLQDLKELKVLVAD KDLLDLHKSLVCTALRGKLVGFSEMEANFKNLSRGLVNVAAKLTHNKDVRDLFVDLVEKVFVEPCRSDHWPLSDVRFLLNQ YSASVHSLDGRHQALWDRYMGLTRGCLLRLYHD

靶点信息

研究背景	Acidic fibroblast growth factor is mitogenic for a variety of different cell types and acts by stimulating mitogenesis or inducing morphological changes and differentiation. The FIBP protein is an intracellular protein that binds selectively to acidic fibroblast growth factor (aFGF). It is postulated that FIBP may be involved in the mitogenic action of aFGF. Two transcript variants encoding different isoforms have been found for this gene.
基因ID	9158
基因名	FIBP
Swiss	O43427
别名	FIBP;FGFIBP;FIBP-1;TROFAS

产品验证



Western blot analysis of FIBP expressed in BxPC-3, Mouse skeletal muscle using FIBP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn