

SND1 Rabbit pAb

货号: AYP19336

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	101kDa
实测分子量	102kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,HL-60,MCF7,Mouse liver,Rat liver
细胞定位	Cytoplasm,Melanosome,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 26-285 of human SND1 (NP_055205.2).
序列	MVLSGCAIIVRGQPRGGPPPERQINLSNIRAGNLARRAAATQPDAKDTDPDEPWAFPAREFLRKKLIGKEVCFTIENKTPQGR EYGMILGKDTNGENIAESLVAEGLATRREGMRANNPEQNRLECEEQAKAAKGMWSEGNCSHTIRDLYTIENPRHFV DSHHQKPVNAIIEHVRDGSVVRALLLPDYLVTVMLSGIKCPTFRREADGSETPEPFAAEAKFFTESRLLQRDVQIILESCHN QNILGTILHPNGNIT

靶点信息

研究背景	This gene encodes a transcriptional co-activator that interacts with the acidic domain of Epstein-Barr virus nuclear antigen 2 (EBNA 2), a transcriptional activator that is required for B-lymphocyte transformation. Other transcription factors that interact with this protein are signal transducers and activators of transcription, STATs. This protein is also thought to be essential for normal cell growth. A similar protein in mammals and other organisms is a component of the RNA-induced silencing complex (RISC). [provided by RefSeq, Jul 2016]
基因ID	27044
基因名	SND1
Swiss	Q7KZF4
别名	SND1;TDRD11;Tudor-SN;p100

产品验证



Western blot analysis of SND1 expressed in Mouse liver, Rat liver, HeLa using SND1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn