

# MLF1 Rabbit pAb

货号: **AYP19365**

## 产品信息

|       |                                                                                                        |
|-------|--------------------------------------------------------------------------------------------------------|
| 反应    | Human,Mouse,Rat                                                                                        |
| 宿主    | Rabbit                                                                                                 |
| 克隆性   | Polyclonal                                                                                             |
| 预测反应  |                                                                                                        |
| 应用    | WB IF/ICC                                                                                              |
| 推荐浓度  | <b>WB:</b> 1:500 - 1:2000<br><b>IF/ICC:</b> 1:50 - 1:200                                               |
| 理论分子量 | 23kDa/27kDa/29kDa/30kDa/33kDa                                                                          |
| 实测分子量 | 38kDa                                                                                                  |
| 形式    | Liquid                                                                                                 |
| 保存条件  | Store at -20°C. Avoid freeze / thaw cycles.<br>Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物   | Unconjugated                                                                                           |
| 阳性对照  | Mouse lung,Rat brain,Rat lung                                                                          |
| 细胞定位  | Cytoplasm,Nucleus                                                                                      |
| 纯化    | Affinity purification                                                                                  |

## 抗原信息


|      |                                                                                                                                                                                                                                                                             |
|------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-268 of human MLF1 ( NP_071888.1).                                                                                                                                                           |
| 序列   | MFRMLNSSFEDDPFFSEILAHRENMRQMIRSFSEPFGRDLLSISDGRGRAHNRGRHNDGEDSLTHTDVSSSQTMDQM VSNMRNYMQKLERNFGQLSVDPNGHSFCSSSVMTYSKIGDEPPKVFQASTQTRRAPGGIKETRKAMRDSDSGLEKMAIG HHIHDRAHVIKSKNKKTGDEEVNQEFINMNESDAHAFDEEWQSEVLKYKPRHNLGNTRMRSVGHENPGSRELKRREK PQQSPAIEHGRRSNVLGDKLHIKGSVKSNNK |

## 靶点信息

|      |                                                                                                                                                                                                                                                                                                                                                  |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 研究背景 | This gene encodes an oncoprotein which is thought to play a role in the phenotypic determination of hemopoietic cells. Translocations between this gene and nucleophosmin have been associated with myelodysplastic syndrome and acute myeloid leukemia. Multiple transcript variants encoding different isoforms have been found for this gene. |
|------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

|       |        |
|-------|--------|
| 基因ID  | 4291   |
| 基因名   | MLF1   |
| Swiss | P58340 |
| 别名    | MLF1   |

## 产品验证



Western blot analysis of MLF1 expressed in Mouse lung,Rat brain,Rat lung using MLF1 Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)