

# HLX Rabbit pAb

货号: **AYP19374**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	50kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60,Mouse lung
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 80-280 of human HLX ( NP_068777.1).
序列	HPHASFQAAARSPLRPTPVVAPSEVPAGFPQRLSPLSAAYHHHPQQQQQQQPQQQPPPPPRAGALQPPASGTRVV PNPHHSGSAPAPSSKDLKFGIDRILSAEFDPKVKEGNTLRDLTSLTGGRPAGVHLSGLQPSAGQFFASLDPINEASAILSP LNSNPRNSVQHQQFQDTFPGPYAVLTKDTMPQTYKRKRWSR

## 靶点信息

研究背景	Transcription factor required for TBX21/T-bet-dependent maturation of Th1 cells as well as maintenance of Th1-specific gene expression. Involved in embryogenesis and hematopoiesis (By similarity).
基因ID	3142

基因名	HLX
Swiss	Q14774
别名	HLX;HB24;HLX1

## 产品验证



Western blot analysis of HLX expressed in HL-60, Mouse lung using HLX Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)