

OSBPL10 Rabbit pAb

货号: **AYP19420**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	77kDa/83kDa
实测分子量	100kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-431,HT-29,A-549,Rat brain
细胞定位	Cytoplasm,cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 200-400 of human OSBP L10 (NP_060254.2).
序列	ASPCSQRHLSVAGPVVTITHHKSPAAARRAKSQYSGQLHEVREMMNQVEGQKNLVHAIESLPGSGPLTALDQDLLLLK ATSAATLSCLGECLNLLQQSVHQAGQPSQKPGASENILGWHGSKSHSTELKNGTLGSLPSASANITWAILPNSAEDEQT SQPEPEPNSGSELVLSSEDEKSDNEDKEETELGVMEQDRSII

靶点信息

研究背景	This gene encodes a member of the oxysterol-binding protein (OSBP) family, a group of intracellular lipid receptors. Like most members, the encoded protein contains an N-terminal pleckstrin homology domain and a highly conserved C-terminal OSBP-like sterol-binding domain. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	114884

基因名	OSBPL10
Swiss	Q9BXB5
别名	OSBPL10;ORP10;OSBP9

产品验证



Western blot analysis of OSBPL10 expressed in Rat brain,A-431,HT-29 using OSBPL10 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane . Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn