

PARVB Rabbit pAb

货号: **AYP19447**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	37kDa/41kDa/45kDa
实测分子量	42kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,A375,HeLa,Mouse heart,Rat heart
细胞定位	Cell junction,Cell membrane,Cell projection,Cytoplasm,Cytoplasmic side,Peripheral membrane protein,Z line,cytoskeleton,focal adhesion,lamellipodium,myofibril,sarcomere
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human PARVB (NP_037459.2).
序列	MSSAPRSPTPRPRRMKKDESFLGKLGGLTARKRRAREVSDLQEEGKNAINSPMSPALVDVHPEDTQLEENEERTMIDPTS KEDPKFKELVKVLLDWINDV

靶点信息

研究背景	This gene encodes a member of the parvin family of actin-binding proteins, which play a role in cytoskeleton organization and cell adhesion. These proteins are associated with focal contacts and contain calponin homology domains that bind to actin filaments. This family member binds to alphaPIX and alpha-actinin, and it can inhibit the activity of integrin-linked kinase. This protein also functions in tumor suppression. Alternative splicing of this gene results in multiple transcript variants.
------	--

基因ID	29780
基因名	PARVB
Swiss	Q9HBI1
别名	PARVB;CGI-56

产品验证



Western blot analysis of PARVB expressed in Mouse heart,Rat heart,U-87MG using PARVB Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn