

# CDK5RAP3 Rabbit pAb

货号: **AYP19459**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	31kDa/46kDa/56kDa/59kDa
实测分子量	65kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG, 293T, LO2
细胞定位	Cytoplasm, Nucleus, centrosome, cytoskeleton, microtubule organizing center
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CDK5RAP3 (NP_788276.1).
序列	MEDHQHVPIDIQTSKLLDWLDRRHCSLKWQSLVLTIREKINAAIQDMPESSEEIAQLLSGSYIHYFHCLRILDLLKGTEASTK NIFGRYSSQRMKDQWEIIALYEKDNTYLVELSSLLVRNVNVEIPSLKKQIAKCCQLQEQEYSRKEEECCAGAAEMREQFYHSC KQYGITGENVRGELLALVKDLPSQLAEIGAAAQQLGAEIDVYQASVGFVCEPTEQVLPMLRFVQKRGNSTVYEWRTGT EPSVVERPHLEELPEQVAEDAIDWGDGVEAVSEGTDGSGISAEAAAGIDWGIFPES

## 靶点信息

研究背景	This gene encodes a protein that has been reported to function in signaling pathways governing transcriptional regulation and cell cycle progression. It may play a role in tumorigenesis and metastasis. A pseudo gene of this gene is located on the long arm of chromosome 20. Alternative splicing results in multiple transcript variants that encode different isoforms.
------	--

基因ID	80279
基因名	CDK5RAP3
Swiss	Q96JB5
别名	CDK5RAP3;C53;HSF-27;IC53;LZAP;MST016;OK/SW-cl.114;PP1553

## 产品验证



Western blot analysis of CDK5RAP3 expressed in U-87MG,293T,LO2 using CDK5RAP3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)