

— ABLYBIO, Help Your Research



VPS25 Rabbit pAb

货号: **AYP19462**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200
理论分子量	20kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,Mouse brain,Mouse heart,Rat liver,Rat lung
细胞定位	Cytoplasm,Endosome membrane,Nucleus,nucleoplasm
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human VPS25 (NP_115729.1).
------	---

靶点信息

研究背景	This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II (ESCRT-II). This protein complex functions in sorting of ubiquitinated membrane proteins during endocytosis. A pseudogene of this gene is present on chromosome 1.
基因ID	84313
基因名	VPS25
Swiss	Q9BRG1 (https://www.uniprot.org/uniprotkb/Q9BRG1/entry)
别名	VPS25,DERP9,EAP20,FAP20,VPS25 Rabbit pAb,Dermal papilla-derived protein 9,ELL-associated protein of 20 kDa,ESCRT-II complex subunit VPS25

产品验证



Western blot analysis of VPS25 expressed in Mouse brain,Rat liver,BxPC-3 using VPS25 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)