

# VPS25 Rabbit pAb

货号: **AYP19462**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB IF/ICC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200
理论分子量	20kDa
实测分子量	20kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BxPC-3,Mouse brain,Mouse heart,Rat liver,Rat lung
细胞定位	Cytoplasm,Endosome membrane,Nucleus,nucleoplasm
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-176 of human VPS25 ( NP_115729.1).
序列	MAMSFEPWPQYRFPPFFTLQPNVDTRQKQLAAWCSLVLSFCRLHKQSSMTVMEAQESPLFNNVKLQRKLPVESIQIVLEE LRKKGNLEWLDKSKSSFLIMWRRPEEWGKLIYQWVSRSGQNNVFTLYELTNGEDTEDEEFHGLDEATLLRALQALQQE HKAEIITVSDGRGVKFF

## 靶点信息

研究背景	This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II ( ESCRT-II). This protein complex functions in sorting of ubiquitinated membrane proteins during endocytosis. A pseudogene of this gene is present on chromosome 1.
基因ID	84313

基因名	VPS25
Swiss	Q9BRG1
别名	VPS25;DERP9;EAP20;FAP20

## 产品验证



Western blot analysis of VPS25 expressed in Mouse brain,Rat liver,BxPC-3 using VPS25 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)