

— ABLYBIO, Help Your Research



# YY1AP1 Rabbit pAb

货号: **AYP19466**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	50kDa/75-97kDa
实测分子量	88kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HT-29,BxPC-3,SGC-7901,HepG2,Mouse liver,Rat lung
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-80 of human YY1AP1 (NP_060723.2).
------	-------------------------------------------------------------------------------------------------------------------

## 靶点信息

研究背景	The encoded gene product presumably interacts with YY1 protein; however, its exact function is not known. Alternative splicing results in multiple transcript variants encoding different isoforms.
------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

基因ID	55249
基因名	YY1AP1
Swiss	Q9H869 ( <a href="https://www.uniprot.org/uniprotkb/Q9H869/entry">https://www.uniprot.org/uniprotkb/Q9H869/entry</a> )
别名	YY1AP1,GRNG,HCCA1,HCCA2,YY1AP,YY1AP1 Rabbit pAb,Hepatocellular carcinoma susceptibility protein,Hepatocellular carcinoma-associated protein 2

## 产品验证



Western blot analysis of YY1AP1 expressed in Mouse liver,Rat lung,HT-29 using YY1AP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)