

# NLN Rabbit pAb

货号: **AYP19471**

## 产品信息

反应	Human,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	80kDa
实测分子量	80kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-549,Rat liver
细胞定位	Cytoplasm,Mitochondrion intermembrane space
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 38-260 of human NLN ( NP_065777.1).
序列	SSYTVAGRNVLRWDLSPQIKTRTEELIVQTKQVYDAVGMLGIEEVTYENCLQALADVEVKYIVERTMLDFPQHVSDDKEV RAASTEADKRLSRFDIEMSMRGDIFERIVHLQETCDLGKIKPEARRYLEKSIKMGKRNGHLHLEPQVQNEIKSMKKRMSELCI DFNKNLNEDDTFLVFSKAELGALPDDFIDSLEKTDKDKYKITLKYPHYFPVMKKCCIPET

## 靶点信息

研究背景	This gene encodes a member of the metallopeptidase M3 protein family that cleaves neurotensin at the P ro10-Tyr11 bond, leading to the formation of neurotensin(1-10) and neurotensin(11-13). The encoded protein is likely involved in the termination of the neurotensinergic signal in the central nervous system and in the gastrointestinal tract.
基因ID	57486

基因名	NLN
Swiss	Q9BYT8
别名	NLN;AGTBP;EP24.16;MEP;MOP

## 产品验证



Western blot analysis of NLN expressed in HeLa, A-549, Rat liver using NLN Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)