

LARP5 Rabbit pAb

货号: **AYP19484**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	80kDa
实测分子量	70KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,A-431,U-251MG
细胞定位	Cytoplasm,cytosol
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 519-738 of human LARP 5 (NP_055970.1).
序列	PPKPPSPSFELGLSSFPLPGAAGNLKTEDLFENRLSSLIIIGPSKERTLSADASVNTLPVVVSREPSVPASCAVSATYERSPSP AHLPPDDPKVAEKQRETHSVDRRLPSALTATAACKSVQVNGAATELRKPSYAEICQRTSKEPPSSPLQPQKEQKPNTVGCGKEE KKLAEPAERYREPPALKSTPGAPRDQRRPAGGRPSPSAMGKRLSREQSTPPKSPQ

靶点信息

研究背景	This gene encodes a member of an evolutionarily conserved protein family implicated in RNA metabolism and translation. Members of this family are characterized by the presence of an La motif, which is often located adjacent to one or more RNA recognition motifs (RRM). Together, the two motifs constitute the functional region of the protein and enable its interaction with the RNA substrate. This protein family is divided into five sub-families: the genuine La proteins and four La-related protein (LARP) sub-families. The protein encoded by this gene belongs to LARP sub-family 4. It is a cytoplasmic protein that may play a stimulatory role in translation.
基因ID	23185
基因名	LARP4B
Swiss	Q92615
别名	LARP4B;KIAA0217;LARP5

产品验证



Western blot analysis of LARP5 expressed in 293T, A-431, U-251MG using LARP5 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn