

# GAB2 Rabbit pAb

货号: **AYP19489**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	69kDa/74kDa
实测分子量	74kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Mouse kidney
细胞定位	Cell membrane,Cytoplasm
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 377-676 of human GAB 2 (NP_536739.1).
序列	ENSRVAATIPRRNTLPAMDNSRLHRASSCETYEYPQRGGESAGRSAESMSDGVGSFLPGKMIVGRSDSTNSEDNYVPM NPGSSTLLAMERAGDNSQSVYIPMSPGAHHFDSLGYPSTTLPVHRGSPSRGSEIQPPPVNRNLKPDRKAKPTPLDLRNNTVI DELPFKSPITKSWSRANHTFNSSSSQYCRPISTQSITSTDSGDSEENYVPMQNPVSASPVPSTNSPAPKKSTGSVDYLAL DFQPSSPSPHRKPSTSSVTSDEKVDYVQVDKEKTQALQNTMQEWTDVQRSSSEPSKGAKL

## 靶点信息

研究背景	This gene is a member of the GRB2-associated binding protein (GAB) gene family. These proteins contain pleckstrin homology (PH) domain, and bind SHP2 tyrosine phosphatase and GRB2 adapter protein. They act as adapters for transmitting various signals in response to stimuli through cytokine and growth factor receptors, and T- and B-cell antigen receptors. The protein encoded by this gene is the principal activator of phosphatidylinositol-3 kinase in response to activation of the high affinity IgE receptor. Two alternatively spliced transcripts encoding different isoforms have been described for this gene.
------	---

基因ID	9846
基因名	GAB2
Swiss	Q9UQC2
别名	GAB2

## 产品验证



Western blot analysis of GAB2 expressed in HeLa, Mouse kidney using GAB2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)