

TLX1 Rabbit pAb

货号: **AYP19563**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:1000 - 1:2000
理论分子量	26kDa/34kDa
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,K562,Mouse spleen,Rat spleen
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human TLX1 (NP_001182446.1).
序列	MEHLGPHHLHPGHAEPISFGIDQILNSPDQGGCMGPASRLQDGEYGLGCLVGGAYTYGGGGSSAAATGAGGAGAYGTG GPGGPGGPAGGGGACSMGPLTGSYVNMALAGGPGPGGGGGSSGGAGALSAAGVIRVPAHRPLAGAVAHPQPLATGL PTVPSVPAMPGVNLTGLTFPVMESNRRYTKDRFTG

靶点信息

研究背景	This gene encodes a nuclear transcription factor that belongs to the NK-linked or NK-like (NKL) subfamily of homeobox genes. The encoded protein is required for normal development of the spleen during embryogenesis. This protein is also involved in specification of neuronal cell fates. Ectopic expression of this gene due to chromosomal translocations is associated with certain T-cell acute lymphoblastic leukemias. Alternative splicing results in multiple transcript variants.
------	---

基因ID	3195
基因名	TLX1
Swiss	P31314
别名	TLX1;HOX11;TCL3

产品验证



Western blot analysis of TLX1 expressed in Mouse spleen, Rat spleen, Jurkat using TLX1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn